VERMEIREN

SLINGS Albatros

INSTRUCTION MANUAL
MODE D'EMPLOI
GEBRUIKSAANWIJZING
GEBRAUCHSANWEISUNG
ISTRUZIONI PER L'USO
MANUAL DE INSTRUCCIONES
INSTRUKCJA OBSŁUGI


Fig 1


Fig 2


Fig 3


Fig 4

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 1.7 Explanation of symbols6

2 Use 7

 2.1 Attaching to the buttons of the patient hoist.....7

 2.2 Raising8

3 Maintenance and inspection 9

 3.1 Regular maintenance and inspection9

 3.2 Shipping and storage9

 3.3 Care.....10

 3.4 Disinfection10

4 Guarantee 13

5 Disposal 13

6 Declaration of conformity 13

7 Maintenance plan 14

8 Disinfection book 14

Preface

First of all we want to thank you for putting your trust in us by selecting one of our products.

The expected lifetime of your sling is strongly influenced by the care and maintenance of the sling.

This manual will help you get acquainted with the operation of your slings. In these document you can also find the maintenance instructions for a longer lifetime of your sling.

Following the user instructions and the maintenance instructions are an essential part of the guarantee.

This manual reflects the latest product developments. The company Vermeiren has the right to introduce changes without the obligation to adapt or replace previously delivered models.

For any further questions, please consult your specialist dealer.

1 Product description

1.1 *Intended use*

⚠ CAUTION: Risk of injury – Make sure the patient can be lifted with slings.

The slings are used for the Vermeiren "Albatros" patient hoist (active patient hoist) and is intended for people with walking difficulties.

The sling is designed to transport 1 person and not to be used for persons younger than 12 years.

The slings allow full use by persons disabled by:

- paresis
- limb deformations
- cachexia
- atrophy
- and also for aged persons.

When providing for individual requirements:

- body size and weight (max. 200kg (440.9 lb.))
- physical and psychological condition
- environment

should be taken into consideration.

The manufacturer is not liable for damage caused by the lack of or improper service or as a result of not following instructions in this manual.

Read the manual of the patient hoist "Albatros" before using the slings. It is very important to understand the content of the manual "Albatros" very well.

Compliance with the user and maintenance instructions are an essential part of the guarantee conditions.

1.2 *Delivery*

The Vermeiren slings shall be delivered with:

- 1 sling
- Manual

Before use check if everything is included and that no products are damaged (example by transport, ...).


1.3 Packaging

The slings are packed separate in a box.

1.4 Drawings


- 1 = Backrest
- 2 = Mounting strap patient hoist
- 3 = Waist belt with velcro
- 4 = Armpit support

1.5 Technical specifications

Choose the right size of sling. The slings are available in two sizes: M, XL.
 Make sure that the chosen model, size and material meets the needs of the patient.

The different sizes are indicated by different colors of the zoom.

- Medium (M) Yellow
- Extra large (XL) Black

Technical specification below are valid for the Vermeiren slings.

Make	Vermeiren	
Address	Vermeirenplein 1/15, B-2920 Kalmthout	
Type	Sling	
Maximum occupant mass	200 kg (440.9 lb.)	
Description		Dimensions
Order sizes	M (yellow)	XL (black)
Weight range person*	40 - 120 kg (88.2 - 264.6 lb.)	100 - 200 kg (220.5 - 440.9 lb.)
Model	Active	
Use	Normal use for active patient	
A	860 mm (33.86 in)	1060 mm (41.73 in)
B	380 mm (14.96 in)	380 mm (14.96 in)
Washing temperature	+ 60 °C (+ 140 °F)	
Storage and use temperature	+ 5 °C (+ 41 °F)	+ 45 °C (+ 113 °F)
Storage and use humidity	30%	70%
<p>We reserve the right to introduce technical changes. Measurement tolerance ± 15 mm / kg (0.59 in. / 3,3 lb.) * The specified patient weight range is a guideline.</p>		

Table 1: Technical Specifications slings


The slings are made of 100% polyester in a solid structure (closed substance). This material is easy to maintain, durable and slides easily.

100% polyester has a limited absorbing and heat insulating capacity so it is a pleasure to wear clothing made of natural fiber.

The slings comply to the requirements set up in:

EN ISO 10535: 2006

1.6 Identification of product label

1.7 Explanation of symbols

Maximum user mass

Machine wash on 60°C /140°F

Do not use chlorine bleach

Do not use a tumble dry

Do not iron

CE conformity

Read the instruction manual before use

Observe the safety instructions

2 Use


This chapter describes the normal use.

The sling and patient hoist may only be used by trained personnel. The care provider is always responsible for the safety of the patient.

2.1 Attaching to the buttons of the patient hoist

The slings shall be attached by the mounting straps. These mounting straps shall be suspended on the buttons of the patient hoist.

The mounting straps of the slings have 3 loops in different colors.

With the 3 different loops it is possible to raise patients from different heights.

Color	Description
Red	Tall patients, the hands of the patient can come close to the grip arms of the patient hoist so the patient can help actively to stand up.
Yellow	Combination
Black	Small patients: Patient is relatively short compared to the height of the buttons.

Table 2: Color codes loops mounting straps

2.2 Raising

Safety instructions

- Inform and train the personnel about the safety and risks when using the patient hoist and sling. This training shall be done by competent persons.
- The responsible person in your organization must decide if one or more care providers are necessary.
- Apply the brakes of the wheelchair, brancard, etc.
- Plan the raising of the patient prior to the actual raising to ensure safe and correct raising.
- Before moving the patient to another place be sure that the patient is safely and comfortably in the sling.
- Don't use damaged loops or slings.
- Check that the loops are attached correctly to the buttons of the patient hoist. Check this before the patient is raised.
- Never leave the patient without supervision.
- Check the maximum working weight of the patient hoist. The smallest working weight between patient hoist and sling is determining.

2.2.1 Raising and lowering instructions

The patient can be raised with the patient hoist by following the instructions below:

1. If the patient is raising from a wheelchair, apply the brakes of the wheelchair and remove the arm supports.
2. Spread out the sling and hold the marking label to you.
3. Place the sling behind the back of the patient. *(Cover page, Fig 1, arrow 1)* Be aware the sling is positioned in the middle.
4. Position the armpit cushions of the sling under the armpits of the patient. *(Cover page, Fig 1)*
5. Mount the waist belt with the velcro strip around the waist of the patient. Do not put on too tight or loose. *(Cover page, Fig 2, arrow 1)*
6. Position the patient lift as close to the chair in the desired position.
7. Position the legs of the patient hoist in the desired position around the chair.
8. Position the footplate of the patient hoist as close to chair. Place the feet of the patient on the footplate.
9. Position the lower legs parallel with the leg support. Adjust the height and depth of the leg support so the leg support will sit comfortable just under the knees.
10. Use the loops to mount the straps to the buttons of the patient hoist. *(Cover page, Fig 3: position 1 can be used for small patients and position 2 for big patients)*
11. Raise the patient by pushing on the manual control.
12. Now the patient can be moved to the desired place with the patient lift. *(Cover page, Fig 4)*

The patient can be lowered with the patient hoist by following the instructions below:

1. If the patient is lowered to a wheelchair, apply the brakes of the wheelchair and remove the arm supports.
2. Position the wheelchair as close to the patient hoist. When lowering a patient into a chair, try to position the patient as close as possible to the back rest of the chair.
3. Lower the patient by pushing on the manual control.
4. Position the patient in the chair.
5. Remove the loops of the mounting straps from the buttons of the patient hoist. *(Cover page, Fig 3)*
6. Remove the waist belt with the velcro strip. *(Cover page, Fig 2, arrow 2)*
7. Remove the sling behind the back of the patient. *(Cover page, Fig 1, arrow 2)*

3 Maintenance and inspection

The expected lifetime of the sling is 2 - 4 years, depending on type of use, storage, regular maintenance, servicing and cleaning/washing.

3.1 Regular maintenance and inspection

The maintenance activities to ensure a good condition for your sling are described below:

- For every use and after washing
 - Check the mounting straps on wear and damage
 - Check the zooms on wear and damage
 - Check the loops on wear and damage
 - Check the backrest on wear and damage

3.2 Shipping and storage

The shipping and storage shall be according following instructions:

- Store in a dry place (between +5 °C / +41°F and +45 °C / +113°F).
- The relative humidity of the air should be between 30% and 70%.
- Components must be shipped or stored without being subjected to damage (not clamping between something, be careful with sharp objects, ...).


3.3 Care

Clean the sling for every new user. Note the instructions below for cleaning the sling:

- Clean the sling with a cloth moistened with hot water.
- Use a mild commercial detergent for removing stubborn dirt.
- Stains can be removed by using a sponge or a soft brush.
- Do not use strong cleaning liquids like solvents, nor use hard brushes.
- Never clean with steam and/or pressure washers.

The sling can also be washed in the washing machine. Below some instructions for washing in the washing machine:

- The slings can be washed at + 60 ° C (+ 140 °F) in the normal machine wash. Do not use chlorine bleach.
- Do not put the slings in the tumble dry.
- Do not iron the slings but just fold them up.

3.4 Disinfection

- ⚠ WARNING: Dangerous Products - The use of disinfectants is restricted to authorized personnel.**
- ⚠ WARNING: Dangerous Products, change on irritate your skin - You should wear suitable protective clothing because the disinfectants could irritate your skin. For this purpose you should also take note of the product information of the solutions concerned.**

All parts of the slings can be treated by scrubbing with a disinfectant.

All steps taken to disinfect rehabilitation equipment, their components or other accessory parts are to be recorded in a disinfection report containing a minimum of the following information (with product documentation appended):

Date of the disinfection	Reason	Specification	Substance and concentration	Signature
--------------------------	--------	---------------	-----------------------------	-----------

Table 3: Example of a disinfection book

Abbreviations used in column 2 (reason):

V = Suspected infection IF = Infection case W = Repetition I = Inspection

For a blank sheet of a disinfection book see § 8.

The recommended disinfectants for scrubbing (based on the list provided by the Robert Koch Institute, RKI) are standing in the table below. The current state of the disinfectants included in the RKI list can be obtained from the Robert Koch Institute (RKI) (homepage: www.rki.de).

Active substance	Product name	Laundry disinfection		Surface disinfection (scrubbing-/wiping disinfection)		Disinfection of excretions 1 part sputum or stools + 2 parts diluted solution or 1 part urine + 1 part diluted solution						Area of effectiveness	Manufacturer or Supplier
		Diluted solution	Time to take effect	Diluted solution	Time to take effect	Sputum		Stools		Urine			
						Diluted solution	Time to take effect	Diluted solution	Time to take effect	Diluted solution	Time to take effect		
		%	Hr.	%	Hr.	%	Hr.	%	Hr.	%	Hr.		
Phenol or phenol derivative	Amocid	1	12	5	6	5	4	5	6	5	2	A	Lysoform
	Gevisol	0,5	12	5	4	5	4	5	6	5	2	A	Schülke & Mayr
	Helipur			6	4	6	4	6	6	6	2	A	B. Braun
	m-cresylic soap solution (DAB 6)	1	12	5	4							A	
	Phenol	1	12	3	2							A	
Chlorine, organic or inorganic substances with active chlorine	Chloramin-T DAB 9	1,5	12	2,5	2	5	4					A ¹ B	
	Clorina	1,5	12	2,5	2	5	4					A ¹ B	Lysoform
	Trichlorol	2	12	3	2	6	4					A ¹ B	Lysoform
Per combinations	Apesin AP100 ²			4	4							AB	Tana PROFESSIONAL
	Dismozon pur ²			4	1							AB	Bode Chemie
	Perform ²			3	4							AB	Schülke & Mayr
	Wofesteril ²			2	4							AB	Kesla Pharma
Formaldehyde and/or other aldehydes or derivatives	Aldasan 2000			4	4							AB	Lysoform
	Antifect FD 10			3	4							AB	Schülke & Mayr
	Antiseptica surface disinfection 7			3	6							AB	Antiseptica
	Apesin AP30			5	4							A	Tana PROFESSIONAL
	Bacillocid special			6	4							AB	Bode Chemie
	Buraton 10F			3	4							AB	Schülke & Mayr
	Desomed A 2000			3	6							AB	Desomed
	Hospital disinfectant cleaner			8	6							AB	Dreiturm
	Desomed Perfekt			7	4							AB*	Desomed
	Formaldehyde-solution (DAB 10), (formaline)	1,5	12	3	4							AB	
	Incidin Perfekt	1	12	3	4							AB	Ecolab
	Incidin Plus			8	6							A	Ecolab
Kohrsolin	2	12	3	4							AB	Bode Chemie	

Active substance	Product name	Laundry disinfection		Surface disinfection (scrubbing-/wiping disinfection)		Disinfection of excretions 1 part sputum or stools + 2 parts diluted solution or 1 part urine + 1 part diluted solution						Area of effectiveness	Manufacturer or Supplier
		Diluted solution	Time to take effect	Diluted solution	Time to take effect	Sputum		Stools		Urine			
						Diluted solution	Time to take effect	Diluted solution	Time to take effect	Diluted solution	Time to take effect		
		%	Hr.	%	Hr.	%	Hr.	%	Hr.	%	Hr.		
Formaldehyde and/or other aldehydes or derivatives	Lysoform	4	12	5	6							AB	Lysoform
	Lysoformin	3	12	5	6							AB	Lysoform
	Lysoformin 2000			4	6							AB	Lysoform
	Melsept	2	12	4	6							AB	B. Braun
	Melsitt	4	12	10	4							AB	B. Braun
	Minutil	2	12	6	4							AB	Ecolab
	Multidor			3	6							AB	Ecolab
	Nüscosept			5	4							AB	Dr. Nüsken Chemie
	Optisept			7	4							AB*	Dr. Schumacher
	Pursept-FD			7	4							AB*	Merz
Ultrasol F	3	12	5	4							AB	Fresenius Kabi	
Amphoteric surfactants (amfotensiden)	Tensodur 103	2	12									A	MFH Marienfelde
Lye	Lime-milk ³							20	6			A ³ B	

1 Not effective against myco-bacteria when service disinfecting, especially in the presence of blood.
2 Not suitable for disinfecting blood-contaminated or porous surfaces (e.g. raw wood).
3 Useless for tuberculosis; preparation of Lime-milk: 1 part dissolved lime (calcium hydroxide) + 3 parts water.
* Checked for effectiveness on viruses in accordance with checking methods of the RKI (Federal Health Reporting 38 (1995) 242).
A: Suitable for killing vegetative bacterial germs including myco-bacteria as well as fungi, including fungal spores.
B: Suitable for deactivating viruses.

Table 4: Disinfectants

Kindly consult your specialist dealer if you have queries on matters related to disinfection; he will gladly assist you.

4 Guarantee

Excerpt from the "General Business Conditions":

(...)

5. The guarantee period for warranty claims is 24 months.

(...)

The guarantee excludes damage arising from structural changes to our products, insufficient maintenance, defective or improper handling or storage or the use of non-original parts.

(...)

5 Disposal

When disposing of the sling, contact your local disposal centre or return the product to your specialist dealer.

6 Declaration of conformity

The manufacturer or his authorized representative :

N.V. VERMEIREN N.V

Address :

Vermeirenplein 1/15

2920 Kalmthout

Belgium

declares under his sole responsibility that the CE marked medical devices :

Productgroup: Slings active patient hoist

Brand: Vermeiren

Type: Active

have been classified as class I, according to annex IX MDD 93/42/EEC, rule 1,

and is manufactured in full conformity with the European instructions below - including the latest modifications - and with the national law, that organizes this directions :

Medical devices directive MDD 93/42/EEC

and is in conformity with the relevant European harmonized standards:

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ENGLISH

CONTRACTUAL WARRANTY

We offer 5 years of warranty on standard wheelchairs, lightweight wheelchairs 4 years. Electronic wheelchairs, tricycles, beds and other products: 2 years (batteries 6 months) and multiposition wheelchairs 3 years. This warranty is limited to the replacement of defective or spare parts.

APPLICATION CONDITIONS

In order to claim this warranty, part "B" of this card has to be given to your official Vermeiren dealer. The warranty is only valid when parts are replaced by Vermeiren in Belgium.

EXCEPTIONS

This warranty is not valid in case of:

- damage due to incorrect usage of the wheelchair,
- damage during transport,
- involvement in an accident,
- a dismount, modification or repair carried outside of our company and/or official Vermeiren dealership,
- normal wear of the wheelchair,
- non-return of the warranty card

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All rights reserved, including translation.
Vermeiren Suisse S.A.
Vermeiren Nederland B.V.