


| Material |                        |         |          |
|----------|------------------------|---------|----------|
| No.      | Plan                   | Article | Quantity |
| 1        | XP02840 (High 100)     | 2994649 | 1        |
|          | XP02841 (High 120)     | 2994650 | 1        |
| 2        | SANMED C-460/5202 Blue | 2994994 | 4,3m     |
| 3*       | Slide fastener         | 1995604 | 4,2m     |
| 4*       | Slide fastener head    | 1995605 | 2        |
| 5*       | Sewing thread          | W200016 | 20m      |




Final dimensions of mattress: 1920x1340x100mm or 1920x1340x120mm /  
Wymiary koncowe materaca: 1920x1340x100mm lub 1920x1340x120mm

EP00469 / A

100 or 120

\* not depicted / nie ujeto na rysunku

| MATTRESS FOR INTERVAL XXL 140 |            |            |            | Plan No.  |
|-------------------------------|------------|------------|------------|---|
|                               |            |            |            | EP00469 / A   |
| Scale: 1:12                   | Drafted    | Verified   | Confirmed  |  |
| Date                          | 29/10/2019 | 04/02/2021 | 04/02/2021 |   |
| Initials                      | AZ         | EB         | TK         |   |