



No.	Plan	Article	Qty
1	FRONT COVER	29036	1
2	DUAL LOCK 3M 15x25mm	2998521	2
3	RUBBER	3992256	4

Before stick 3M Dual Lock prepare surface on plastic part:

1. Roughen the surface by abrasive pad
2. Clean surface with Isopropyl Alcohol (IPA), propan-2-OL.
3. Press 3M Dual Lock after stick

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

SV SYNC FRONT COVER SUB-ASSY			
Scale: 1:5	Drafted	Verified	Confirmed
Date	29/08/2025	14/10/2025	14/10/2025
Initials	BK	BK	MK

Plan No. E00108126/A.2