



#	Reference	Art number	Description	Option	Qty
1	EC_WIRE6	2998659	EC WIRE EON/SYNC MAIN POWER CABLE		1
2	EC_L5, EC_L8	1920896	EC L SV KUARZ/SYNC/EON REAR/BRAKE LIGHT		2
3	EC_L4, EC_L7	1920895	EC L SV KUARZ/SYNC/EON REAR BLINKER		2
4	EC_COM10	3992118	SV CARPO2 25K RESISTOR ADJUSTABLE BH		1
4	EC_L1, EC_L3	1920894	EC L SV KUARZ/SYNC/EON FRONT BLINKER		2
5	EC_BAT1, EC_BAT2	1917939	EC BAT BS E BATTERIJ 12V DAB12-44EV-HD AGM		2
6	EC_WIRE7	2998661	EC WIRE SV EON/SYNC CABLE FOR REAR LIGHTS		1
7	EC_WIRE5	2998660	EC WIRE SV EON/SYNC CONTROL CABLE		1
8	EC_WIRE4	2998658	EC WIRE SV EON MAIN POWER FUSE 120A		1
9	EC_WIRE2	1921858	EC WIRE KUARZ CABLE FOR FRONT LIGHTS		1
10	EC_SW1	3992376	EC SW BATTERY SWITCH 150A A701		1
11	EC_MOT1	3992041	EC MOT CERES SE MOTORSET 450W	10 Km/h	1
		1922754	EC MOT SV SYNC MOTORSET 500W	12 Km/h	1
		3993577	EC MOT SV CERES SE MOTORSET 650W	15 Km/h	1
12	EC_L2, EC_L9	1920897	EC L SV SYNC/EON FRONT LIGHT		2
13	EC_COM6	3992597	EC COM EON/SYNC PCB		1
14	EC_COM9	3990622	EC COM SV CONTROLLER S-DRIVE 70	10 Km/h 12 Km/h	1
		3990623	SV CARPO2 CONTROLLER S-DRIVE 120	15 Km/h	1
15	EC_COM8	1921863	EC COM SV KUARZ CHARGE HOLDER		1
16	EC_COM3	3992014	EC COM SV CARPO2 BUZZER BH		1
17	EC_COM2	3992019	EC COM TAC05 5K VR + ALU BLOCK BH I		1
18	EC_COM1	1922287	SV KUARZ KEY SETS BH 2801		1
19	UNDEPICTED	2998581	CABLE TIES AHM PM-150		15
20	EC_WIRE9	2998662	EC WIRE SV FRONT LIGHT SPLITTER		1

BK BLACK
 RD RED
 LBL LIGHT BLUE
 BL BLUE
 LGR LIGHT GREEN
 GR GREEN
 VI VILLET
 OR ORANGE
 YE YELLOW
 WH WHITE
 LRD LIGHT RED
 PK PINK
 GY GRAY
 / AND
 I OR

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
 PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
 WHOLE OR IN PART, OR USED FOR ANY REASON OF MANUFACTURE,
 EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORIZATION
 FROM VERMEIREN GROUP

BOM EC SYNC

Size: A4	Drafted	Verified	Conformed
Date:	2026-04-14	2026-04-14	2026-04-14
Initials:	FM	PM	PM

2998655/A

Plan No.

2998655/A

