



#	Reference	Art number	Description	Optrion Qty
1	EC_WIRE6	2998659	EC WIRE EON/SYNC MAIN POWER CABLE	1
2	EC_L5, EC_L8	1920896	EC L SV KUARZ/SYNC/EON REAR/BRAKE LIGHT	2
3	EC_L4, EC_L7	1920895	EC L SV KUARZ/SYNC/EON REAR BLINKER	2
4	EC_L1, EC_L3	1920894	EC L SV KUARZ/SYNC/EON FRONT BLINKER	2
5	EC_BAT1, EC_BAT2	1910262	EC BAT BS E BATTERIJ 38AH	38AH 2
		1904684	EC BAT BS E BATTERIJ 38AH	50AH 2
6	EC_WIRE7	2998661	EC WIRE SV EON/SYNC CABLE FOR REAR LIGHTS	1
7	EC_WIRE5	2998660	EC WIRE SV EON/SYNC CONTROL CABLE	1
8	EC_WIRE4	2998658	EC WIRE SV EON MAIN POWER FUSE 120A	1
9	EC_WIRE2	1921858	EC WIRE KUARZ CABLE FOR FRONT LIGHTS	1
10	EC_SW1	3992376	EC SW BATTERY SWITCH 150A A701	1
11	EC_MOT1	3992043	EC MOT SV ERIS MOTORSET 320W	1
12	EC_L2	1920897	EC L SV SYNC/EON FRONT LIGHT	1
13	EC_COM6	3992597	EC COM EON/SYNC PCB	1
14	EC_COM9	3990622	EC COM SV CONTROLLER S-DRIVE 70	1
15	EC_COM8	1921863	EC COM SV KUARZ CHARGE HOLDER	1
16	EC_COM3	3992014	EC COM SV CARPO2 BUZZER BH	1
17	EC_COM2	3992019	EC COM TACOS 5K VR + ALU BLOCK BH I	1
18	EC_COM1	1922287	SV KUARZ KEY SETS BH 2801	1
19	UNDEPICTED	2998581	CABLE TIES AHM PM-150	5
20	EC_COM10	3992118	SV CARPO2 25K RESISTOR ADJUSTABLE BH	1

- BK BLACK
- RD RED
- LBL LIGHT BLUE
- BL BLUE
- LGR LIGHT GREEN
- GR GREEN
- VI VILLET
- OR ORANGE
- YE YELLOW
- WH WHITE
- LRD LIGHT RED
- PK PINK
- GY GRAY
- / AND
- I OR

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY PURPOSE OF MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORIZATION FROM VERMEIREN GROUP

BOM EC EON

Size: A4 Drafted Verified Conformed

Date: 2025-10-14

Initials: FM

2998657/A

Plan No. 2998657/A

