


No.	Plan	Article	Qty
1	P00103622	2998665	1
2	P00098666	2998629	1
3	P00100623	29136	1
4	WHEEL LINE PROFILE	1922057	1
5	1922805	1922805	2
6	3991137	3991137	1
7	WASHER M8 DIN127B	1001570	3
8	3991119	3991119	1
9	NUT DIN985 M8	3991743	2
10	NUT M6 DIN985 BLACK	3991110	4
11	SCREW M6X16 ISO7380 BLACK	2994067	4
12	WASHER M6 BLACK	2994068	2
13*	WHEEL BLOCK PROFILE	1922054	1

* - Optional equipment, non depicted, block profile wheel

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORIZATION FROM VERMEIREN GROUP.

EON RIGHT FRONT WHEEL				Plan No. E00085893/A.1
Scale: 1:3	Drafted	Verified	Confirmed	
Date	22/07/2025	11/10/2025	11/10/2025	
Initials	BK	BK	MK	