


Material			
No.	Plan	Article	Quantity
1	XP02603 (High 100)	2994388	1
	XP02839 (High 120)	2994648	1
2	SANMED C-460/7000 Grey (Vinterno)	2994993	4,3m
	SANMED C-460/5202 Blue (Vermeiren)	2994994	4,3m
3*	Slide fastener (Vinterno)	2994224	4,2m
	Slide fastener (Vermeiren)	1995604	4,2m
4*	Slide fastener head (Vinterno)	2994225	2
	Slide fastener head (Vermeiren)	1995605	2
5*	Sewing thread (Vinterno)	W200319	20m
	Sewing thread (Vermeiren)	W200016	20m




Final dimensions of mattress: 1920x1165x100mm or 1920x1165x120mm /
Wymiary końcowe materaca: 1920x1165x100mm lub 1920x1165x120mm

EP00419 / C

100 or 120

* not depicted / nie ujeto na rysunku

MATTRESS FOR INTERVAL XXL 120 WITH CENTRAL BRAKE / INTERVAL 120				Plan No.
				EP00419 / C
Scale: 1:10	Drafted	Verified	Confirmed	
Date	30/07/2019	04/02/2021	04/02/2021	
Initials	AZ	EB	TK	