




ECLIPS X4/KIDS SIDE FRAME WITH B05 2005				Plan No.
				E000213/I
Scale: 6:25	Drafted	Verified	Confirmed	
Date	24/03/2005	18/12/2024	18/12/2024	
Init.	SD	PM	AS	

Material							
No.	Plan	Article	Number	No.	Plan	Article	Number
1	Y005296	15503..	1	31	X003852	1670082	1
	Y005322	15403..	1	32	X000404	G000789	2
2	Y005293	15502..	1	33	M6x45 DIN912/8.8	1003109	5
3	Y005047 (L)	15332..	1	34	---	1005860	1
	Y005048 (R)	15333..	1	35	X000850	1244904	1
4	Y005989 (L)	15635..	1	36	---	1001602	1
	Y005990 (R)	15636..	1	37	X000861	1245004	1
5	Y005553	1551154	1	38	---	1001680	1
6	---	1914611	1	39	---	1992600	1
7	X002636	1994438	1	40	---	1001561	1
8	X002670	1547400	1	41	---	1994605	1
9	X003384	1997149	1	42	X002132	1992423	1
10	---	1992559	2	43	---	1992597	1
11	X003303 (L)	1596382	1	44	X002169	1992791	1
	X003302 (R)	1596282	1	45	X002163	1488100	1
12	X002205	1489300	1	46	X002130 (L)	1992421	1
13	Y006648	1595182	1		X002131 (R)	1992422	1
14	X002749	1553400	1	47	Y005129	G084204	2
15	---	1001684	1	48	M16 type S DIN17222	1911653	1
16	---	1006850	3	49	---	1006110	1
17	---	1001610	1	51	---	1001682	1
18	X002135	1992405	1	52	X003480	1996670	1
19	---	1001564	13	53	---	1001570	2
20	Y006196	1595004	1	54	M6x35 DIN912/8.8	1001496	4
21	X003481	1996587	2	55	---	1905244	1
22	---	1008492	13	56	XC00037	3990106	1
23	X001674	1990601	1	-	-	-	1
24	X002679	1994452	1				
25	---	2990070	2				
26	X000582	1005740	1				
27	---	1004550	1				
28	---	1001589	1				
29	M6x35	1006057	1				
30	X002437	1992699	1				

ECLIPS X4/KIDS SIDE FRAME WITH B05 2005				Plan No.
				E000213/I
Scale: 6:25	Drafted	Verified	Confirmed	
Date	24/03/2005	18/12/2024	18/12/2024	
Init.	SD	PM	AS	