



NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

SV CERES SE MAIN BASIC			
Scale:	1:5	Drafted	Verified
Date	14/05/2016	07/04/2023	07/04/2023
Initials	KZ	PM	TK

3998132

Plan No.
EC00324/K

No.	Plan	Article	Qty
1	YC04059	1C99204	1
2A*	YC02054	1C74704	1
2B*	P00057276	2674904	1
3	YC01726	1C56304	1
4	YC01723	1C56004	1
5	YC01722	1C55904	1
6	YC01721	1C55800	1
7	YC01720DP	1C55704	1
8	YC01701	1C55504	2
9	XC00716	3993272	1
10	XC00715	1C15404	1
11	XC00538	3992261	1
12	XC00452	GC00097	2
13	XC00274	3990732	2
14	XC00179	3990713	4
15	X004436	1903928	2
16	X001016	1001682	4
17	P00051750	2672204	1
18	P00051694	2672304	1
19	Washer	3992610	2
20	---	3992314	1
21	M8x90	3992232	2
22	Spring	3991958	2
23	Nut M14	3991097	1
24	M6x25	3991095	2
25	---	3990781	1

26	M6x45	1991017	4
27	12V 38Ah	1909446	2
28	M6x16	1008459	4
29	M8x35	1007516	2
30	M6x20	1006999	6
31	Screw D3.9x16	1005767	2
32	M10x55	1005636	1
33	Nut M10	1001683	1
34	Nut M6	1001680	16
35	Nut M6	1001670	1
36	Washer	1001564	2
37	Washer	1001554	2

2A* & 2B* - depends of batteries type

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

3998132

SV CERES SE MAIN BASIC

Plan No.

EC00324/K

Scale:	1:5	Drafted	Verified	Confirmed
Date	14/05/2016	07/04/2023	07/04/2023	
Initials	KZ	PM	TK	

