



Bill of Material			
No.	Plan	Article	Quantity
1	---	0990781 (L)	1
	---	0990782 (R)	1
	---	0990783 (L)	1
	---	0990784 (R)	1
2	Y002672	1436704	2
3	Y002673	1436804	2
4	---	F002039	2
5	M6x16 ISO7380	1008459	8
6	M6 DIN985	1001680	4
7	M6x35 DIN933	1001433	2
8	M6x25 DIN963	1001479	2
9	M6 DIN125A	1001564	4
10	YP00948	1835854 (1)	1
	YP00949	1835954 (2)	1
	YP00951	1836054 (3)	1
	YP00952	1836154 (4)	1
11	X003698	1004550	2
12	M6x25 ISO7380	1006998	2

- (1) B03 Right
- (2) B03 Left
- (3) B02/B05 Right |
- (4) B02/B05 Left |

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

B112 ARMREST FOR HEMIPLEGIA

Scale:	2:5	Drafted	Verified	Confirmed
Date	9/10/2012	24/09/2025	24/09/2025	
Initials	NV	AZ	AS	

Plan No.
E000842/C
