



Material			
No.	Plan	Article	Quantity
1	Y017126	1740054	1
2	Ø22	1908818	2
3	Ø3x10	2992305	4
4	25x25	1001617	2
5	VGf40 M8x45	2990016	1
6	Y023056	2051554	1
7	VGf40 M8x16	2990024	1
8	20x300	2993635	1
9	Ø5x16 DIN7337A	1901989	1

E000649 / M

2991004

B31 CRUTCH HOLDER				Plan Number
				E000649/M
Scale: 1:4	Drafted	Verified	Confirmed	
Date	13/09/2010	31/05/2021	31/05/2021	
Initials	TT	TT	PJ	