



BOM			
No.	Plan	Article	Qty
1	Y018809	1776254	1
2	Y017126	1740054	1
3	Y023056	2051554	1
4	-	-	-
5	25x25	1001617	2
6	D22	1908818	2
7	D4,8X12	1001585	1
8	D3x10	1908849	4
9	M8x35 DIN912	2994318	1
10	M8x45	2990016	1
11	M8x16	2990024	1
12	---	2993635	1
13	-	-	-

Adaptor kit 1919334 = 1 + 9

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

1919334 /
1919504

Plan No.
E00049016/F.1

B31 CRUTCH HOLDER TYPE ONE

Scale:	1:4	Drafted	Verified	Confirmed
Date	8/12/2022	27/04/2026	27/04/2026	
Initials	SD	KA	AZ	