



No.	Plan	Article	Qty
1	Basket	3990644	1
2	3990791	3990791	1
3	YC00616	1C31904	1
4	YC00868	1C31804	1
5	Screw M5x16 ISO7380	1006112	4
6	3992273	3992273	2
7	P00125349	GP33900	2
8	Screw M6x50 DIN7991	1001486	2
9	Screw M5X30 DIN912	1903983	2
10	P00049955	T001093	1

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

SV SYNC/EON FRONT BASKET				---
				Plan No.
				E00125348/A.1
Scale:	1:5	Drafted	Verified	Confirmed
Date	21/04/2026	21/04/2026	21/04/2026	21/04/2026
Initials	BK	MK	MK	MK

