



NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Sheet 1/2

2993391

### Forest 3+ Rear Lower Frame

Plan No.

E001058/D

|          |            |            |            |           |
|----------|------------|------------|------------|-----------|
| Scale:   | 1:5        | Drafted    | Verified   | Confirmed |
| Date     | 23.11.2016 | 11.08.2025 | 11.08.2025 |           |
| Initials | SD         | PM         | MK         |           |



| Material |            |         |          |     |               |         |          |
|----------|------------|---------|----------|-----|---------------|---------|----------|
| No.      | Plan       | Article | Quantity | No. | Plan          | Article | Quantity |
| 1        | X001016    | 1001682 | 4        | 18  | M10 DIN985    | 1001683 | 2        |
| 2        | X002282    | 1992854 | 4        | 19  | M8x45 DIN912  | 1003110 | 4        |
| 3        | X002792    | 1555954 | 4        | 20  | M6x12 ISO7380 | 1007000 | 9        |
| 4        | X003403    | 1617400 | 4        | 21  | M5x35 DIN912  | 1007961 | 2        |
| 5        | X004778    | 1908255 | 4        | 22  | M12 DIN125    | 1993362 | 2        |
| 6        | X005576    | 1912735 | 1        | 23  | ---           | 1906791 | 2        |
| 7        | X005609    | 1912734 | 1        | 24  | ---           | 1906792 | 2        |
| 8        | X005637    | 1912761 | 1        | 25  | ---           | 1908165 | 2        |
| 9        | X005638    | 1912762 | 1        | 26  | M5x10         | 1909204 | 2        |
| 10       | Y020969 DP | 20363.. | 1        | 27  | 80A           | 1909646 | 2        |
| 11       | Y021565    | 2025854 | 1        | 28  | 150A          | 1910488 | 1        |
| 12       | Y022051    | 2025954 | 1        | 29* | ---           | 1912309 | 1        |
| 13       | M6 DIN125  | 1000973 | 6        | 30  | LINX 120A     | 1912415 | 1        |
| 15       | M5 DIN125  | 1001550 | 2        | 31  | M6            | 1912736 | 3        |
| 16       | D28        | 1001606 | 2        | 32  | M10x65 DIN912 | 1912777 | 2        |
| 17       | M5 DIN985  | 1001679 | 2        | 33  | ---           | 1999843 | 2        |
|          |            |         |          | 34  | M5x10 ISO7380 | 2991862 | 4        |

\* Option SE50

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

Sheet 2/2

2993391

Forest 3+ Rear Lower Frame

Plan No.

E001058/D

|          |            |            |            |           |
|----------|------------|------------|------------|-----------|
| Scale:   | 1:5        | Drafted    | Verified   | Confirmed |
| Date     | 23.11.2016 | 11.08.2025 | 11.08.2025 |           |
| Initials | SD         | PM         | MK         |           |

