



No.	Plan	Article	Qty
1	P00123153	29234	1
2	TUBE PLUG 30x30 Black	1003198	6
3	SCREW M8X45 DIN912	2995397	4
4	Star Knob VGF 40 M8x20	1001580	1

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

SV SYNC OPTIONAL EQUIPMENT BAR			
Scale:	1:3	Drafted	Verified
Date	26/03/2026	26/03/2026	26/03/2026
Initials	MK	PM	MK

-
Plan No. E00124025/A.1
